


## TPO-11 Dehydration of 2-methylcyclohexanol



Place 20 mL of 2-methylcyclohexanol (**1**) and a stirring rod into a 50 mL Claisen-Vigreux flask topped by an additional funnel in which 5 mL of 85% phosphoric acid (handle with care) are introduced. Warm up gently and stir.

Add the phosphoric acid dropwise. Distil out the products by heating the mixture gently. *Do not allow the temperature of the distillate to rise above 115 °C.*

Transfer the distillate to a separating funnel. Wash the organic layer with 10 mL of cold water, then with 10 mL of 5% sodium hydrogencarbonate and then with as many water necessary.

Transfer the organic layer to a clean, dry erlenmeyer flask and dry it.

Purify your product with the appropriate method.

Analyse your product mixture by GC (gas chromatography) and IR (infrared spectroscopy).

### Cautions

- 85% phosphoric acid is corrosive to skin and clothing. Neutralise any spills with sodium bicarbonate and quickly wash the affected areas of skin with lots of soap and water.
- Methylcyclohexenes are flammable, volatile and have a noxious odor ; be sure to handle them under your student hood. Keep all containers of methylcyclohexene capped.

### Physical data



- 2-methylcyclohexanol (**1**) : FW = 114.2 g.mol<sup>-1</sup> ; d = 0.93 ; bp 165-168°C ; purity : 99% ; very sol. in ethanol.
- 1-methylcyclohexene (**2**) : FW = 96.2 g.mol<sup>-1</sup> ; d = 0.81 ; bp 110°C ; insoluble in water, soluble in diethylether.
- 3-methylcyclohexene (**3**) : d = 0.80 ; bp 104 °C ;
- 4-methylcyclohexene (**4**) : d = 0.79 ; bp 102 °C ;
- methylenecyclohexane (**5**) : d = 0.80 ; bp 102 °C. All insol. in water, very sol. in diethylether.




**Results**


1. Calculate the yield of your synthesis R.
2. Analyse your chromatogram (qualitatively, and then quantitatively with the normalization method).
3. Calculate the yield  $R_c$  corrected by your chromatography analysis for each product obtained.
4. Interpret your experimental IR spectra.
5. IR spectra of 2-methylcyclohexanol are recorded as : spectrum **1** a liquid film and spectrum **2** as  $\text{CCl}_4$  solution. Interpret the difference between these two spectra.
6. Assign the peaks of the  $^1\text{H}$  NMR and  $^{13}\text{C}$  NMR spectra **3** and **4** below to specific groups of protons and carbon atoms on the molecule.
7. Attribute the MS spectrum **5**, and NMR spectra **6** and **7**.

**Questions**

1. Does the E2 mechanism explain the formation of compounds (4) and (5) ?
2. What geometric requirement of the E2 mechanism causes the cis isomer to dehydrate faster ? If the trans isomer undergoes E2 elimination, which product is formed ?
3. How does Saitzev's Rule apply to this experiment ?
4. Show that E1 mechanism accounts for the formation of the by-products (4) and (5).
5. Explain the role of phosphoric acid in this reaction.
6. Explain what would happen if we used hydrochloric acid instead of phosphoric acid ?
7. Explain why the temperature measured at the top of the column (first distillation) must not exceed  $115\text{ }^\circ\text{C}$  ?
8. Explain the difficulty to purify correctly the 1-methylcyclohexene.

**Spectrum 1 : FT-IR Spectrum of (1)****Spectrum 2 : FT-IR Spectrum of (1)**

**Spectrum 3 : Interpret the Carbon NMR Spectrum of (2)****Spectrum 4 : Interpret the Proton NMR Spectrum of (2)****Spectrum 5 : Mass Spectrum of Unknown Compound**

**Spectrum 6 : Carbon NMR Spectrum of (3), (4) or (5) ?****Spectrum 7: Proton NMR Spectrum of (3), (4) or (5) ?**